


Anti-Peroxiredoxin 3 (Prdx3) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P30048

Item Number: CM0005461



Introduction

Rabbit monoclonal antibody against human Peroxiredoxin 3 (PRDX3), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse PRDX3; mitochondrial antioxidant protein.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 3 (PRDX3) Rabbit mAb
Remarks/Alias	AOP1; MER5; AOP-1; PPPCD; SP-22; HBC189; SCAR32; PRO1748; prx-III; Peroxiredoxin 3 (PRDX3)
Species	Human
GeneID (Human)	10935
GeneID	10935
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 157-256 of human Peroxiredoxin 3 (PRDX3) (P30048)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30048
Protein Weight	28kDa
Shipping	Ice bag