


Anti-Peroxiredoxin 3 (Prdx3) Rabbit Polyclonal Antibody For Wb, Elisa - P30048

Item Number: CM0020412



Introduction

CamelBio offers a rabbit anti-PRDX3 antibody validated for WB and ELISA. It detects human and rat PRDX3, a mitochondrial antioxidant enzyme critical in oxidative stress and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 3 (PRDX3) Rabbit pAb
Remarks/Alias	AOP1; MER5; AOP-1; PPPCD; SP-22; HBC189; SCAR32; PRO1748; prx-III
Species	Human
Gene ID	10935
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 63-256 of human Peroxiredoxin 3 (PRDX3) (NP_006784.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P30048
Protein Weight	28 kDa
Shipping	Ice bag