


Anti-Peroxiredoxin 4 (Prdx4) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13162

Item Number: CM0011375



Introduction

Rabbit monoclonal antibody against human Peroxiredoxin 4 (PRDX4), validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PRDX4 (~31 kDa). Ideal for oxidative stress and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 4 (PRDX4) Rabbit mAb
Remarks/Alias	PRX-4; AOE372; AOE37-2; HEL-S-97n; Peroxiredoxin 4 (PRDX4)
Species	Human
Gene ID	10549
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Peroxiredoxin 4 (PRDX4) (Q13162)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13162
Protein Weight	31kDa
Shipping	Ice bag