


Anti-Peroxiredoxin 4 (Prdx4) Rabbit Polyclonal Antibody For Wb, Elisa - Q13162

Item Number: CM0014179



Introduction

Rabbit polyclonal antibody targeting human Peroxiredoxin 4 (PRDX4), validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for oxidative stress and antioxidant enzyme research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Peroxiredoxin 4 (PRDX4) Rabbit pAb
Remarks/Alias	PRX-4; AOE372; AOE37-2; HEL-S-97n; 4)
Species	Human
GeneID	10549
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-271 of human Peroxiredoxin 4 (PRDX4) (NP_006397.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13162
Protein Weight	31kDa
Shipping	Ice bag