


Anti-Peroxiredoxin 5 (Prdx5) Rabbit Monoclonal Antibody For Wb, Elisa - P30044

Item Number: CM0009859



Introduction

Rabbit monoclonal antibody against human PRDX5, validated for WB and ELISA. Detects a ~22kDa band in human, mouse, and rat samples. Targets the mitochondrial antioxidant enzyme peroxiredoxin-5, protecting cells from oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 5 (PRDX5) Rabbit mAb
Remarks/Alias	PLP; ACR1; B166; PRXV; PMP20; PRDX6; prx-V; SBB110; AOEB166; HEL-S-55; Peroxiredoxin 5 (PRDX5)
Species	Human
GeneID (Human)	25824
GeneID	25824
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30044
Protein Weight	22kDa
Shipping	Ice bag