


Anti-Peroxiredoxin 5 (Prdx5) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P30044

Item Number: CM0006905



Introduction

Rabbit monoclonal antibody against human Peroxiredoxin 5 (PRDX5) validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for oxidative stress research. High specificity targeting mitochondrial, cytoplasmic, and peroxisomal isoforms.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Peroxiredoxin 5 (PRDX5) Rabbit mAb
Remarks/Alias	PLP; ACR1; B166; PRXV; PMP20; PRDX6; prx-V; SBBI10; AOEB166; HEL-S-55; Peroxiredoxin 5 (PRDX5)
Species	Human
GeneID	25824
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 53-214 of human Peroxiredoxin 5 (PRDX5)(NP_036226.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30044
Protein Weight	22kDa
Shipping	Ice bag