


Anti-Peroxiredoxin 5 (Prdx5) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30044

Item Number: CM0029988



Introduction

High-quality rabbit polyclonal antibody against Peroxiredoxin 5 (PRDX5), validated for Western blot, IF/ICC, and ELISA. Suitable for studying oxidative stress and antioxidant defense in human, mouse, and rat samples.

[Learn More](#)

Parameter	Specification
Product Name	Peroxiredoxin 5 (PRDX5) Rabbit pAb
Remarks/Alias	PLP; ACR1; B166; PRXV; PMP20; PRDX6; prx-V; SBBI10; AOEB166; HEL-S-55; Peroxiredoxin 5 (PRDX5)
Species	Human
GeneID (Human)	25824
GeneID	25824
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 53-214 of human Peroxiredoxin 5 (PRDX5) (NP_036226.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30044
Protein Weight	22kDa
Shipping	Ice bag