

Anti-Peroxiredoxin 6 Monoclonal Antibody For Wb, Elisa - P30041

Item Number: CM0005085

Introduction

Rabbit monoclonal antibody against human Peroxiredoxin 6 (PRDX6), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for oxidative stress and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 6 Rabbit mAb
Remarks/Alias	PRX; p29; AOP2; 1-Cys; NSGPx; aiPLA2; LPCAT-5; HEL-S-128m; Peroxiredoxin 6
Species	Human
Gene ID (Human)	9588
Gene ID	9588
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30041
Protein Weight	25kDa
Shipping	Ice bag