

Anti-Peroxiredoxin 6 Polyclonal Antibody For Wb - P30041

Item Number: CM0019622



Introduction

Rabbit polyclonal antibody against human Peroxiredoxin 6 (PRDX6), for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat. 25 kDa protein. Antioxidant enzyme with peroxidase and phospholipase activity.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 6 Rabbit pAb
Remarks/Alias	PRX; p29; AOP2; 1-Cys; NSGPx; aiPLA2; LPCAT-5; HEL-S-128m; Peroxiredoxin 6
Species	Human
Gene ID (Human)	9588
Gene ID	9588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human Peroxiredoxin 6 (NP_004896.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30041
Protein Weight	25kDa
Shipping	Ice bag