


Anti-Peroxiredoxin 6 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Elisa - P30041

Item Number: CM0005195



Introduction

This rabbit monoclonal antibody against human Peroxiredoxin 6 (KD validated) detects a ~25 kDa protein, suitable for WB, IHC-P, and ELISA with human and rat cross-reactivity. Ideal for oxidative stress and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Peroxiredoxin 6 Rabbit mAb
Remarks/Alias	PRX; p29; AOP2; 1-Cys; NSGPx; aiPLA2; LPCAT-5; HEL-S-128m
Species	Human
GeneID	9588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human Peroxiredoxin 6 (NP_004896.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P30041
Protein Weight	25 kDa
Shipping	Ice bag