

Anti-Perp Rabbit Monoclonal Antibody For Wb, Ihc, If - Q96Fx8

Item Number: CM0005944

Introduction

PERP Rabbit Monoclonal Antibody (Q96FX8) validated for WB, IHC, IF, ELISA in human, mouse, rat. Targets p53 apoptosis effector related to PMP-22, involved in desmosome junctions and epithelial integrity. Ideal for apoptosis and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	PERP Rabbit mAb
Remarks/Alias	THW, KCP1, EKVP7, OLMS2, PIGPC1, KRTCAP1, dj496H19.1, PERP
Species	Human
GeneID (Human)	64065
GeneID	64065
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 94-193 of human PERP (Q96FX8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FX8
Protein Weight	21kDa
Shipping	Ice bag