

Anti-Pex11B Rabbit Polyclonal Antibody For Wb And Elisa - O96011

Item Number: CM0019426

Introduction

Rabbit polyclonal antibody against human PEX11B (Peroxin-11B), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studies on peroxisomal proliferation and division.

[Learn More](#)

Parameter	Value
Product Name	PEX11B Rabbit pAb
Remarks/Alias	PEX14B; PEX11beta; PEX11-BETA; PEX11B
Species	Human
Gene ID	8799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human PEX11B (NP_003837.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96011
Protein Weight	28kDa
Shipping	Ice bag