

Anti-Pex12 Rabbit Polyclonal Antibody For Western Blot & Elisa - O00623

Item Number: CM0017342

Introduction

This PEX12 rabbit polyclonal antibody is validated for Western blot and ELISA, targeting human peroxisome assembly protein 12 (PAF-3). Ideal for studies of peroxisome biogenesis and PEX5 receptor recycling. Order from CamelBio for reliable IVD raw materials.

[Learn More](#)

Parameter	Value
Product Name	PEX12 Rabbit pAb
Remarks/Alias	PAF-3; PBD3A; PEX12
Species	Human
Gene ID (Human)	5193
Gene ID	5193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-359 of human PEX12 (NP_000277.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00623
Protein Weight	41kDa
Shipping	Ice bag