

Anti-Pex14 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075381

Item Number: CM0029702



Introduction

Rabbit polyclonal antibody targeting human PEX14, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Essential for studies on peroxisomal protein import and peroxisome biogenesis disorders.

[Learn More](#)

Parameter	Value
Product Name	PEX14 Rabbit pAb
Remarks/Alias	NAPP2; PBD13A; Pex14p; dj734G22.2; PEX14
Species	Human
Gene ID (Human)	5195
Gene ID	5195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-377 of human PEX14 (NP_004556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075381
Protein Weight	41kDa
Shipping	Ice bag