


Anti-Pex19 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P40855

Item Number: CM0028415



Introduction

Rabbit polyclonal antibody against human PEX19 (Peroxin-19) validated for WB, IF-P, ELISA. Detects PEX19 in human, mouse, rat. Ideal for peroxisomal biogenesis and membrane protein trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	PEX19 Rabbit pAb
Remarks/Alias	PXF; HK33; PMP1; PMPI; PXMP1; PBD12A; D1S2223E; PEX19
Species	Human
GeneID (Human)	5824
GeneID	5824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human PEX19 (NP_002848.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40855
Protein Weight	33kDa
Shipping	Ice bag