

Anti-Pex2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P28328

Item Number: CM0020559

Introduction

PEX2 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects PEX2 in human, mouse, rat. Essential for peroxisome biogenesis and PEX5 recycling research. Suitable for pexophagy and peroxisomal disorder studies.

[Learn More](#)

Parameter	Specification
Product Name	PEX2 Rabbit pAb
Remarks/Alias	PAF1; PMP3; ZWS3; PBD5A; PBD5B; PMP35; PXMP3; RNF72; PEX2
Species	Human
GeneID (human)	5828
GeneID	5828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human PEX2 (NP_001165558.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28328
Protein Weight	35kDa
Shipping	Ice bag