


Anti-Pex3 Polyclonal Antibody For Wb, Elisa - P56589

Item Number: CM0029262



Introduction

Rabbit polyclonal antibody against human PEX3 (TRG18/PBD10A). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for peroxisome biogenesis and Zellweger syndrome research.

[Learn More](#)

Parameter	Value
Product Name	PEX3 Rabbit pAb
Remarks/Alias	TRG18; PBD10A; PBD10B; PEX3
Species	Human
GeneID (human)	8504
GeneID	8504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 144-373 of human PEX3 (NP_003621.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56589
Protein Weight	42kDa
Shipping	Ice bag