

Anti-Pex5 Polyclonal Antibody For Wb, If/Icc, Elisa - P50542

Item Number: CM0018103



Introduction

Rabbit polyclonal antibody against human PEX5, validated for WB, IF/ICC, and ELISA. Cross-reactive with mouse and rat. Targets peroxisomal targeting signal 1 receptor; key for peroxisome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PEX5 Rabbit pAb
Remarks/Alias	PXR1; PBD2A; PBD2B; PTS1R; RCDP5; PTS1-BP; PEX5
Species	Human
GeneID (Human)	5830
GeneID	5830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 364-631 of human PEX5 (NP_000310.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50542
Protein Weight	71kDa
Shipping	Ice bag