

Anti-Pex6 Rabbit Pab For Wb, Elisa - Q13608

Item Number: CM0013538

Introduction

Rabbit polyclonal antibody targeting human PEX6, a peroxisomal AAA ATPase essential for peroxisome protein import. Validated in WB and ELISA; cross-reacts with mouse and rat. Ideal for peroxisome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PEX6 Rabbit pAb
Remarks/Alias	PAF2; HMLR2; PAF-2; PBD4A; PDB4B; PXAAA1; PEX6
Species	Human
Gene ID (Human)	5190
Gene ID	5190
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 741-980 of human PEX6 (NP_000278.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13608
Protein Weight	104kDa
Shipping	Ice bag