

Anti-Pfas Polyclonal Antibody For Wb, Elisa - O15067

Item Number: CM0019505

Introduction

Rabbit polyclonal antibody targeting human PFAS (PURL, FGAMS) for Western blot and ELISA; cross-reacts with mouse and rat; suitable for purine biosynthesis research and assay development.

[Learn More](#)

Parameter	Value
Product Name	PFAS Rabbit pAb
Remarks/Alias	PURL; FGAMS; GATD8; FGARAT; FGAR-AT; PFAS
Species	Human
GeneID (human)	5198
GeneID	5198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PFAS (NP_036525.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15067
Protein Weight	145kDa
Shipping	Ice bag