


Anti-Pfdn1 Polyclonal Antibody For Wb, Elisa - O60925

Item Number: CM0025798



Introduction

Rabbit polyclonal antibody against PFDN1, a prefoldin subunit involved in protein folding. Validated for WB and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	PFDN1 Rabbit pAb
Remarks/Alias	PDF; PFD1; PFDN1
Species	Human
GeneID (Human)	5201
GeneID	5201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human PFDN1 (NP_002613.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60925
Protein Weight	14kDa
Shipping	Ice bag