


Anti-Pfdn2 Polyclonal Antibody For Wb, Elisa - Q9Uhv9

Item Number: CM0012991



Introduction

Rabbit polyclonal antibody against human PFDN2 (Prefoldin subunit 2). Validated for WB and ELISA, cross-reacts with human and mouse. A useful tool for chaperone and protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	PFDN2 Rabbit pAb
Remarks/Alias	PFD2; PFDN2
Species	Human
GeneID (human)	5202
GeneID	5202
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHV9
Protein Weight	17kDa
Shipping	Ice bag