


Anti-Pfdn5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99471

Item Number: CM0029874



Introduction

Anti-PFDN5 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat PFDN5, a prefoldin chaperone regulating MYC transcription. Ideal for protein folding and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PFDN5 Rabbit pAb
Remarks/Alias	MM1; MM-1; PFD5; PFDN5
Species	Human
GeneID (Human)	5204
GeneID	5204
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99471
Protein Weight	17kDa
Shipping	Ice bag