


Anti-Pfdn5 Rabbit Monoclonal Antibody For Wb, If, Elisa - Q99471

Item Number: CM0008385



Introduction

High-specificity rabbit monoclonal antibody against human PFDN5 (MM1), validated for WB, IF, ELISA. Cross-reacts with mouse and rat. Ideal for studying c-Myc modulation and chaperone functions.

[Learn More](#)

Parameter	Specification
Product Name	PFDN5 Rabbit mAb
Remarks/Alias	MM1; MM-1; PFD5; PFDN5
Species	Human
GeneID (human)	5204
GeneID	5204
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-154 of human PFDN5 (Q99471).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99471
Protein Weight	17 kDa
Shipping	Ice bag