


Anti-Pfkfb3 Rabbit Monoclonal Antibody (Ipfk-2) For Wb, Ihc-P, Elisa - Q16875

Item Number: CM0007052



Introduction

Rabbit monoclonal antibody against human PFKFB3 (iPFK-2), validated for Western Blot, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for glycolysis and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PFKFB3 Rabbit mAb
Remarks/Alias	PFK2; IPFK2; iPFK-2; PFKFB3
Species	Human
Gene ID (Human)	5209
Gene ID	5209
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human PFKFB3 (NP_004557.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16875
Protein Weight	57kDa/58kDa/59kDa/60kDa
Shipping	Ice bag