


Anti-Pfkfb3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16875

Item Number: CM0017371



Introduction

High-quality rabbit polyclonal antibody targeting PFKFB3 (6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 3). Validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Ideal for metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	PFKFB3 Rabbit pAb
Remarks/Alias	PFK2; IPFK2; iPFK-2; PFKFB3
Species	Human
GeneID (Human)	5209
GeneID	5209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PFKFB3 (NP_004557.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16875
Protein Weight	60kDa
Shipping	Ice bag