


Anti-Pfkfb3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16875

Item Number: CM0030611



Introduction

High-quality Anti-PFKFB3 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, and rat. Suitable for research on fructose-2,6-bisphosphate metabolism. - Q16875

[Learn More](#)

Parameter	Value
Product Name	PFKFB3 Rabbit pAb
Remarks/Alias	PFK2; IPFK2; iPFK-2; PFKFB3
Species	Human
Gene ID (Human)	5209
Gene ID	5209
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16875
Protein Weight	60kDa
Shipping	Ice bag