


Anti-Pfkl Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P17858

Item Number: CM0007250



Introduction

High-quality Anti-PFKL Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets liver-type phosphofructokinase, a key glycolysis enzyme. Suitable for metabolism and immunology research.

[Learn More](#)

Parameter	Value
Product Name	PFKL Rabbit mAb
Remarks/Alias	PFK-B, PFK-L, ATP-PFK, PFKL
Species	Human
GeneID	5211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 680-780 of human PFKL (NP_002617.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17858
Protein Weight	85kDa
Shipping	Ice bag