

Anti-Pfkm Polyclonal Antibody For Wb, If/Icc, Elisa - P08237

Item Number: CM0027290



Introduction

Rabbit polyclonal antibody against human PFKM (ATP-dependent 6-phosphofructokinase, muscle type) validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen for glycolysis research.

[Learn More](#)

Parameter	Value
Product Name	PFKM Rabbit pAb
Remarks/Alias	GSD7; PFK1; PFKA; PFKX; PFK-1; PFK-A; ATP-PFK; PPP1R122; PFKM
Species	Human
GeneID (Human)	5213
GeneID	5213
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 501-780 of human PFKM (NP_001160160.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08237
Protein Weight	85kDa
Shipping	Ice bag