

Anti-Pfkm Rabbit Monoclonal Antibody For Wb - P08237

Item Number: CM0012703



Introduction

Rabbit monoclonal antibody targeting human ATP-dependent 6-phosphofructokinase muscle type (PFKM/P08237). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for glycolysis and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	PFKM Rabbit mAb
Remarks/Alias	GSD7; PFK1; PFKA; PFKX; PFK-1; PFK-A; ATP-PFK; PPP1R122; PFKM
Species	Human
GeneID (human)	5213
GeneID	5213
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Fructose 6 Phosphate Kinase (P08237).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08237
Protein Weight	85k Da
Shipping	Ice bag