


Anti-Pfcp Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q01813

Item Number: CM0025500



Introduction

High-quality rabbit polyclonal antibody against PFKP, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat ~86 kDa PFKP. Ideal for glycolysis and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PFKP Rabbit pAb
Remarks/Alias	PFKF; PFK-C; PFK-P; ATP-PFK; PFKP
Species	Human
Gene ID (Human)	5214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 336-540 of human PFKP (NP_002618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01813
Protein Weight	86kDa
Shipping	Ice bag