


Anti-Pfcp Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q01813

Item Number: CM0003810



Introduction

Rabbit monoclonal antibody targeting human PFKP, validated for Western blot, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for glycolysis and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PFKP Rabbit mAb
Remarks/Alias	PFKF; PFK-C; PFK-P; ATP-PFK; PFKP
Species	Human
Gene ID	5214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 336-540 of human PFKP (NP_002618.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01813
Protein Weight	86kDa
Shipping	Ice bag