


Anti-Pfcp Rabbit Polyclonal Antibody For Wb And Elisa - Q01813

Item Number: CM0020683



Introduction

Rabbit polyclonal antibody targeting human PFKP (phosphofructokinase, platelet type). Suitable for WB and ELISA applications. Useful for studying glycolysis, cancer metabolism, and related signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PFKP Rabbit pAb
Remarks/Alias	PFKF; PFK-C; PFK-P; ATP-PFK; PFKP
Species	Human
Gene ID (Human)	5214
Gene ID	5214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 605-784 of human PFKP (NP_002618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q01813
Protein Weight	86kDa
Shipping	Ice bag