

Anti-Pfn2 Rabbit Polyclonal Antibody For Wb, Elisa - P35080

Item Number: CM0027530



Introduction

Rabbit polyclonal antibody targeting human Profilin-2 (PFN2), validated for Western blot (WB) and ELISA. Detects endogenous PFN2 in human, mouse, and rat samples. Ideal for studying actin dynamics and cytoskeleton organization.

[Learn More](#)

Parameter	Value
Product Name	PFN2 Rabbit pAb
Remarks/Alias	PFL; D3S1319E; PFN2
Species	Human
Gene ID (Human)	5217
Gene ID	5217
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-140 of human PFN2 (P35080)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35080
Protein Weight	15kDa
Shipping	Ice bag