

Anti-Pga3 Rabbit Polyclonal Antibody For Wb, Elisa - P0Djd8

Item Number: CM0022800

Introduction

Validated for WB and ELISA, this rabbit polyclonal antibody targets human PGA3 (Pepsinogen-3) with cross-reactivity to mouse and rat. Ideal for digestive enzyme and protease research.

[Learn More](#)

Parameter	Value
Product Name	PGA3 Rabbit pAb
Remarks/Alias	PGA3
Species	Human
GeneID (Human)	643834
GeneID	643834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-180 of human PGA3 (NP_001073275.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DJD8
Protein Weight	42kDa
Shipping	Ice bag