

Anti-Pga5 Rabbit Polyclonal Antibody For Wb, Elisa - P0Djd9

Item Number: CM0018155

Introduction

Anti-PGA5 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat PGA5. Ideal for studying pepsinogen-5 function in digestion and proteolysis.

[Learn More](#)

Parameter	Value
Product Name	PGA5 Rabbit pAb
Remarks/Alias	Pg5; PGA5
Species	Human
GeneID (Human)	5222
GeneID	5222
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-180 of human PGA5 (NP_055039.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DJD9
Protein Weight	42kDa
Shipping	Ice bag