


Anti-Pgam1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P18669

Item Number: CM0011207



Introduction

High-quality rabbit monoclonal antibody targeting human PGAM1, suitable for Western blot, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat. An optimal choice for glycolysis research, cancer metabolism studies, and diagnostic raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	PGAM1 Rabbit mAb
Remarks/Alias	PGAMA; PGAM-B; HEL-S-35; PGAM1
Species	Human
GeneID (Human)	5223
GeneID	5223
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18669
Protein Weight	29kDa
Shipping	Ice bag