

Anti-Pgam2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P15259

Item Number: CM0028333

Introduction

Rabbit polyclonal antibody against human PGAM2, a muscle-specific phosphoglycerate mutase involved in glycolysis. Validated for WB, IF-P, ELISA applications with cross-reactivity to human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	PGAM2 Rabbit pAb
Remarks/Alias	GSD10; PGAMM; PGAM-M; PGAM2
Species	Human
GeneID (human)	5224
GeneID	5224
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human PGAM2 (NP_000281.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15259
Protein Weight	29kDa
Shipping	Ice bag