

Anti-Pgam5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Hs1

Item Number: CM0015463



Introduction

Rabbit polyclonal antibody against human PGAM5, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, and rat. Suitable for studying mitochondrial dynamics, mitophagy, and necroptosis.

[Learn More](#)

Parameter	Value
Product Name	PGAM5 Rabbit pAb
Remarks/Alias	BXLBV68; PGAM5
Species	Human
Gene ID (Human)	192111
Gene ID	192111
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-146 of human PGAM5 (NP_001164014.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96HS1
Protein Weight	32kDa
Shipping	Ice bag