

Anti-Pgam5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Hs1

Item Number: CM0006283



Introduction

Rabbit monoclonal PGAM5 antibody (Bcl-XL-binding protein v68) for WB, IHC-P, IF/ICC, ELISA. Detects human PGAM5 (Q96HS1), a mitochondrial serine/threonine phosphatase regulating mitophagy and necrosis.

[Learn More](#)

Parameter	Value
Product Name	PGAM5 Rabbit mAb
Remarks/Alias	BXLBV68; PGAM5
Species	Human
GeneID	192111
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-146 of human PGAM5 (NP_001164014.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96HS1
Protein Weight	32kDa/28KDa
Shipping	Ice bag