

Anti-Pgc1 Alpha Rabbit Polyclonal Antibody - 070343

Item Number: CM0020827



Introduction

Rabbit polyclonal antibody targeting PGC1 alpha (PPARGC1A), validated for WB, IF, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studying metabolic regulation and mitochondrial biology.

[Learn More](#)

Parameter	Value
Product Name	PGC1 alpha Rabbit pAb
Remarks/Alias	Pgc1; PGC-1; Pgco1; Gm11133; Ppargc1; Pgc-1alpha; A830037N07Rik; PPARGC-1-alpha; PGC1 alpha
Species	Mouse
Gene ID (Human)	10891
Gene ID	19017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse PGC1 alpha (NP_032930.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	070343
Protein Weight	91 kDa
Shipping	Ice bag