


Anti-Pgc1 Beta Monoclonal Antibody For Wb And Elisa - Q86Yn6

Item Number: CM0010768



Introduction

Rabbit recombinant monoclonal anti-PGC1 beta antibody validated for WB and ELISA, reactive with human and rat. Suitable for studying transcriptional coactivator functions in metabolism and mitochondrial biogenesis.

[Learn More](#)

Parameter	Value
Product Name	PGC1 beta Rabbit mAb
Remarks/Alias	PERC; ERRL1; PGC1B; PGC-1(beta); PGC1 beta
Species	Human
GeneID (Human)	133522
GeneID	133522
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 924-1023 of human PGC1 beta (NP_573570.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q86YN6
Protein Weight	113kDa
Shipping	Ice bag