


Anti-Pgc1A Monoclonal Antibody For Wb & Elisa - Q9UbK2

Item Number: CM0003656



Introduction

Recombinant rabbit monoclonal antibody targeting human PGC1 α , validated for Western blot and ELISA. Ideal for studying metabolic regulation, mitochondrial biogenesis, and circadian rhythms. Cross-reacts with human samples.

[Learn More](#)

Parameter	Specification
Product Name	PGC1 α Rabbit mAb
Remarks/Alias	LEM6; PGC1; PGC1A; PGC-1v; PPARGC1; PGC-1alpha; PGC-1(alpha); PGC1 α
Species	Human
GeneID (Human)	10891
GeneID	10891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-480 of human PGC1 α (NP_037393.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBK2
Protein Weight	91kDa
Shipping	Ice bag