


## Anti-Pgc1A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9UbK2

Item Number: CM0014513



### Introduction

Rabbit polyclonal antibody targeting PGC1α (PPARGC1A), suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. For metabolism, gluconeogenesis, and circadian rhythm research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PGC1α Rabbit pAb                                                                                                   |
| Remarks/Alias    | LEM6; PGC1; PGC1A; PGC-1v; PPARGC1; PGC-1alpha; PGC-1(alpha); PGC1α                                                |
| Species          | Human                                                                                                              |
| Gene ID          | 10891                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 610-798 of human PGC1α (NP_037393.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q9UBK2                                                                                                             |
| Protein Weight   | 91kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |