


Anti-Pgc1A/B Rabbit Monoclonal Antibody (Mab) For Wb, Elisa - Q9UbK2/Q86Yn6

Item Number: CM0002716



Introduction

Rabbit monoclonal antibody targeting PGC1α/β (PPARGC1A/PPARGC1B), a transcriptional coactivator regulating metabolism and mitochondrial function. Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for metabolic, thermogenesis, and circadian research.

[Learn More](#)

Parameter	Value
Product Name	PGC1α/β Rabbit mAb
Remarks/Alias	LEM6; PGC-1(alpha); PGC-1alpha; PGC-1v; PGC1; PGC1A; PPARGC1; ERRLL1; PERC; PGC-1(beta); PGC1B; PGC1α/β
Species	Human
GeneID (human)	10891/133522
GeneID	10891/133522
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 699-798 of human PGC1α/β (NP_037393.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBK2/Q86YN6
Protein Weight	113kDa
Shipping	Ice bag