

Anti-Pgd Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P52209

Item Number: CM0003371



Introduction

Anti-PGD rabbit monoclonal antibody validated for WB, IF/ICC, IP, ELISA applications. Detects human, mouse, rat PGD (53 kDa). Essential for pentose phosphate pathway and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PGD Rabbit mAb
Remarks/Alias	6PGD; PGD
Species	Human
GeneID (Human)	5226
GeneID	5226
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 384-483 of human PGD (P52209).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52209
Protein Weight	53kDa
Shipping	Ice bag