

Anti-Pgd Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52209

Item Number: CM0026518



Introduction

Rabbit polyclonal antibody targeting human PGD (6-phosphogluconate dehydrogenase). Validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous PGD in human, mouse, and rat samples. Essential for pentose phosphate pathway research and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PGD Rabbit pAb
Remarks/Alias	6PGD; PGD
Species	Human
GeneID (Human)	5226
GeneID	5226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 184-483 of human PGD (NP_002622.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52209
Protein Weight	53kDa
Shipping	Ice bag