

Anti-Pge Receptor Ep1 (Ptger1) Rabbit Polyclonal Antibody For Wb, Elisa - P34995

Item Number: CM0016995



Introduction

Rabbit polyclonal antibody targeting human PTGER1 (EP1) validated for WB and ELISA. Cross-reacts with mouse. Suitable for studies on renal function, smooth muscle contraction, and GPCR signaling.

[Learn More](#)

Parameter	Value
Product Name	PGE receptor EP1 (PTGER1) Rabbit pAb
Remarks/Alias	EP1; PGE receptor EP1 (PTGER1)
Species	Human
GeneID (Human)	5731
GeneID	5731
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PGE receptor EP1 (PGE receptor EP1 (PTGER1)) (NP_000946.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P34995
Protein Weight	42kDa
Shipping	Ice bag