


Anti-Pge Receptor Ep3 (Ptger3) Rabbit Polyclonal Antibody For Wb, Elisa - P43115

Item Number: CM0027744



Introduction

Rabbit polyclonal antibody against human PTGER3 (EP3 receptor), validated for WB and ELISA with cross-reactivity to mouse and rat. Explore prostaglandin signaling research with CamelBio's reliable raw materials.

[Learn More](#)

Parameter	Value
Product Name	PGE receptor EP3 (PTGER3) Rabbit pAb
Remarks/Alias	EP3; EP3e; EP3-I; EP3-II; EP3-IV; EP3-VI; PGE2-R; EP3-III; Inc003875; PGE receptor EP3 (PTGER3)
Species	Human
GeneID (Human)	5733
GeneID	5733
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-53 of human PGE receptor EP3 (PTGER3) (NP_942009.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43115
Protein Weight	43kDa
Shipping	Ice bag