

Anti-Pge Receptor Ep4 (Ptger4) Rabbit Pab For Wb, Elisa - P35408

Item Number: CM0020959



Introduction

Rabbit polyclonal antibody targeting prostaglandin E2 receptor EP4 subtype (PTGER4). Validated for Western blot and ELISA, with rat cross-reactivity. Useful for studying GPCR-mediated smooth muscle relaxation and renal function.

[Learn More](#)

Parameter	Value
Product Name	PGE receptor EP4 (PTGER4) Rabbit pAb
Remarks/Alias	EP4; EP4R; PGE receptor EP4 (PTGER4)
Species	Human
GeneID (Human)	5734
GeneID	5734
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-488 of human PGE receptor EP4 (PTGER4) (P35408)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P35408
Protein Weight	53kDa
Shipping	Ice bag